




Our Aim To Provide Quality Instruments



PROFESSIONAL DIAGNOSTIC INSTRUMENTS



VI-20-00



VI-20-03



VI-20-04



VI-20-01



VI-20-02



Newest LED technology
One bulb is sufficient

Blade Size:

- No. 1
- No. 2
- No. 3
- No. 4



LED BULB

F.O. Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Integrated Fiber Optics
- MEGALIGHT
- \varnothing 3-4 mm, 5500 Single Fibbers Minimum
- Can Be Treated In Autoclave
- 2.5, 3.5 V Battery Power
- Require 2 x C Size Batteries

FLAPLIGHT FIBER OPTIC LARYNGOSCOPE



Newest LED technology
One bulb is sufficient

Blade Size:

- Mac. 3
- Mac. 4



VI-21-00



LED BULB

F.O. Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Replacable Light Guide
- \varnothing 3-4 mm, 5500 Single Fibbers Minimum
- Can Be Treated In Autoclave
- 2.5, 3.5 V Battery Power
- Require 2 x C Size Batteries
- Swivelling blade tip facilitates difficult incubation.
- The Flaplight is a Macintosh laryngoscope blade with a tip which is adjustable through 70%



VI-22-00



VI-22-03



VI-22-04



VI-22-01



VI-22-02



Newest LED technology
One bulb is sufficient

Blade Size:

- No. 00
- No. 0
- No. 1
- No. 2
- No. 3
- No. 4



LED BULB


F.O. Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Integrated Fiber Optics
- MACINTOSH
- \varnothing 3-4 mm, 5500 Single Fibbers Minimum
- Can Be Treated In Autoclave
- 2.5, 3.5 V Battery Power
- Require 2 x C Size Batteries

MILLER FIBER OPTIC LARYNGSCOPE



VI-23-00



VI-23-05



VI-23-04



VI-23-02



VI-23-03



VI-23-01



Newest LED technology
One bulb is sufficient

Blade Size:

- No. 00
- No. 0
- No. 1
- No. 2
- No. 3
- No. 4



LED BULB

F.O. Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Integrated Fiber Optics
- MILLER
- \varnothing 3-4 mm, 5500 Single Fibbers Minimum
- Can Be Treated In Autoclave
- 2.5, 3.5 V Battery Power
- Require: 2 x C Size Batteries


VI-24-00

VI-24-04

VI-24-05

VI-24-02

VI-24-01

VI-24-03

LED BULB
Conventional Blades:


- Mat Polished Stainless Steel
- Prevents Back Reflections
- Can Be Sterilized
- Compatible According To DIN ISO 7376
- Halogen and LED Bulb 2.5 V
- Require: 2 x C Size Batteries

Blade Size:

- No. 00
- No. 0
- No. 1
- No. 2
- No. 3
- No. 4

MILLER CONVENTIONAL LARYNGOSCOPE

VI-25-00

VI-25-03

VI-25-04

VI-25-01

VI-25-02
Blade Size:

- No. 00
- No. 0
- No. 1
- No. 2
- No. 3
- No. 4


LED BULB
Conventional Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Can Be Sterilized
- Compatible According To DIN ISO 7376
- Halogen and LED Bulb 2.5 V
- Require: 2 x C Size Batteries



VI-26-00



VI-26-03



VI-26-01



VI-26-04



VI-26-02

Blade Size:

- No. 1
- No. 2
- No. 3
- No. 4

Flange with an inter molded durable plastic locking mechanism.
The light source utilizes a crystal-clear plastic tubular fiber optic bundle style designed to reflect bright white light from the handle to the distal tip of the blade.

Economy Conventional Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Can Be Sterilized
- Compatible According To DIN ISO 7376
- Halogen and LED Bulb 2.5 V
- Require: 2 x C Size Batteries

MILLER ECONOMY CONVENTIONAL LARYNGOSCOPE



VI-27-00



VI-27-03



VI-27-01



VI-27-04



VI-27-02

Blade Size:

- No. 1
- No. 2
- No. 3
- No. 4

Flange with an inter molded durable plastic locking mechanism.
The light source utilizes a crystal-clear plastic tubular fiber optic bundle style designed to reflect bright white light from the handle to the distal tip of the blade.

Economy Miller Blades:

- Mat Polished Stainless Steel
- Prevents Back Reflections
- Can Be Sterilized
- Compatible According To DIN ISO 7376
- Halogen and LED Bulb 2.5 V
- Require: 2 x C Size Batteries


VI-28-01

Attachment clip with integrated on/off switch. Secure. Switches off automatically when returned to the pocket.


VI-28-02

Conventional
-Standard illumination
-With 2.5 V vacuum lamp

Fiber Optic
-High quality fiber optics
-High power 2.5 V xenon
-2.5 V LED lamp

-10 single use mini speculum (2.5 and 4mm) with storage bottle packing in black pouch
-Require: 2x AA batteries


FIBER OPTIC ADULT OTOSCOPES

VI-29-01


-Bayonet fitting
-Illumination light 2.5v LED bulb
-Metal handle C
-Plastic Head CPP
-Rheostat dimmer
-05 x Ø 2.5 mm
-05 x Ø 4.0 mm
-Require: 2 x C batteries


VI-29-02

Parker F.O. Otoscope set high quality professional Metal head, chrome plated battery handle Otoscope with extra bright light 2.5v LED bulb illumination provided thought fiber optic system. High quality big swiveling viewing window with 3x magnification lens out put for pneumatic membrane test.

-10 single use mini speculum (Ø 2.5 and 4mm) with storage bottle packing in black pouch.
-Require: 2 x C batteries



VI-30-01

- Bayonet connecting system head fitting ideal for pneumatic tympanic membrane test, 2.5v LED white light
- 3 speculum (size Ø 2.5 3.5 and 4.5mm)
- 3x High quality magnifying lens
- Available: Polish and Matt finish
- Require: 2 x C batteries



VI-30-02

- 2.5v standard illumination LED bulb. Chrome plated handle and head 4x magnifying glasses.
- Reusable speculum size Ø 2.5 3.5 and 4mm include. Connecting for pneumatic Otoscope. Battery handle with rheostat for Adjustment of illumination.
- 3 reusable ear specula.
- Available: Polish and Matt finish.
- Require: 2 x C batteries.

DERMATOSCOPES



VI-31-01

- Supplies: Soft Black pouch packing with extra 0 - 10 scale glass plate
- Require: 2 x C size batteries.
- Available: Polish and Matt finish



VI-31-02

- With approx. 10 times magnification
- focus adjustment – 6 up to +3.5 dpt
- halogen bulb that can be dimmed
- High capacity 2.5 V halogen lamp
- large contact disc (25 mm) with scale setting



VI-32-01 VI-32-02 VI-32-03

Fiber Optic Battery Handles



VI-33-01 VI-33-02 VI-33-03

Conventional Battery Handles

CONVENTIONAL LAYRNGOSCOPE SETS



VI-34-01
Macintosh Conventional
1 Handle + 3 Blades



VI-34-02
Miller Conventional
1 Handle + 4 Blades



VI-35-01
Miller Fiber Optic
2 Handles + 5 Blades



VI-35-02
Macintosh Fiber Optic
1 Handle + 3 Blades

ENT DIAGNOSTIC SET



Set Boxes



VI-36-01
ENT Set 2.5 V
With Accessories

WeChat



VAKMAN INDUSTRIES



Jinnah Islamia Collage Road
Sialkot - 51310 Pakistan.



+92 52 4593350

+92 52 4594276

+92 52 4593849



+92 52 4588276



info@vakman.com.pk



www.vakman.com.pk